

IISER/PUR/1237/AL/SC/25-26

16.12.2025

CORRIGENDUM

Sub: Supply, installation and commissioning of Modula Photoluminescence Spectrometer : reg

Ref: Tender ID: 2025_IISRT_885860_1

1. Based on the clarification request from bidder the revised technical specification is placed at Annexure 1 to this corrigendum.
2. The due dates are extended as follows:-

Bid submission end date: 22.12. 2025 (1500 Hrs)
Bid opening date: 23.12.2025 (1530 Hrs)
3. All other Terms and Conditions remain the same. Bidders may quote accordingly

Thanking You,

Yours Faithfully

Anura Sadath
16/12/25

Assistant Registrar (P&S)



Annexure 1 to corrigendum No IISER/PUR/1237/AL/SC/25-26

dated 16 Dec 2025

Technical Specifications:

S. No	Modular Photoluminescence Spectrometer with following capabilities: Steady-state and time-resolved photoluminescence measurements Fluorescence, Phosphorescence, Time-correlated single-photon counting (TCSPC) (Excitation, emission, constant wavelength synchronous, and constant energy synchronous spectral scanning; 3D excitation/emission scans)	
	Earlier Specs	Modified Specs
1	<p>Sources: ps-Diode Lasers with preferable optical output of pulse width of <100 ps, power >100 mW, and repetition rate up to 100 MHz.</p> <p>Excitation wavelengths: 300 ± 10 nm, 375 ± 10 nm, 405 ± 10 nm, 440 ± 10 nm, 470 ± 10 nm, 510 ± 10 nm, 590 ± 10 nm</p> <p>Auto selection of repetition rate to suit selected time range.</p>	<p>Sources: ps-Diode Lasers with preferable optical output of pulse widths <100 ps, and repetition rate: 2.5 kHz to 20 MHz.</p> <p>Excitation wavelengths: 375 ± 10 nm, (peak power ~120 mW at 10 MHz) 405 ± 10 nm, (peak power ~100 mW at 10 MHz) 440 ± 10 nm, (peak power ~50 mW at 10 MHz) 470 ± 10 nm, (peak power ~80 mW at 10 MHz) 510 ± 10 nm, (peak power ~80 mW at 10 MHz) 635 ± 10 nm, (peak power ~30 mW at 10 MHz)</p> <p>280 ± 10 nm (pulse width ~1ns, avg power ~4 microW)</p> <p>Auto/Manual selection of repetition rate to suit selected time range.</p>
2	<p>Monochromator: Czerny-Turner design with plane gratings for optimized focus at all wavelengths and minimum stray light (double monochromator for excitation and emission) Focal length: ~320 mm Bandpass: 0-25 nm, continuously variable</p>	<p>Monochromator: Czerny-Turner design with plane gratings for optimized focus at all wavelengths and minimum stray light (double monochromator for excitation and emission) Focal length: ~320 mm Bandpass: 0.01-13 nm, continuously variable</p>
3	<p>Detectors: Detection Range: 200 nm to ~1700 nm. Photomultiplier detector(s) with cooled housing (TE or LN₂) and stabilized. Spectral bandwidth: 1.0 to 20 nm or better (continuously variable and computer controlled) Signal to noise ratio: >20,000:1</p>	<p>Detectors: Detection Range: 200 nm to ~1700 nm. Photomultiplier detector(s) with cooled housing (TE or LN₂) and stabilized. Spectral bandwidth: 0.01 to 13 nm or better (continuously variable and computer controlled) Signal to noise ratio: >20,000:1</p>



Annexure 1
16/12/25

	Dark count: <80 cps	Dark count: <80 cps for visible and <80000 cps for NIR
4	Sample Chamber: Helium Cryostat (closed loop)	Sample Chamber: I. Helium Cryostat (closed loop) Temperature range: 10-400 K, for solid samples, Vertically freely standing with lifting mechanism to raise and lower samples when other measurements are done, Ultra-high vacuum pumping station (10^{-7} torr) II. Cuvette pairs: 3.5 mL- 2 Nos III. Solid sample holder/accessory
5	Software/Measurements: <ul style="list-style-type: none"> Time-resolved emission spectra (TRES) for Fluorescence and Phosphorescence Lifetime kinetics Time-resolved anisotropy Excitation-Emission Maps One package for system control, data collection and analysis Automated polarizers for anisotropy measurements (Steady State & TCSPC) Suitable Peltier temperature controller for measuring PL of liquid samples in the range of 0 deg C to 100 deg C. Steady-state Absorption 	Software/Measurements: No change
6	Suitable Computer: Laptop or PC with monitor, keyboard, and mouse for the operations and experiments. Software should be pre-installed, user friendly, original licensed Windows based operating Backup software in the form of CD ROM/pendrive, it should have built in features like facility for scan application, kinetics.	Suitable Computer: No change
7	Installation and Training:	Installation and Training: No change

Amar Sadath
16/12/25



Ph.-0471 2778019,
EMAIL: purchasestores@iisertvm.ac.in

MARUTHAMALA.P.O VITHURA. P.O
THIRUVANANTHAPURAM 695551,
KERALA, INDIA
GST No.32AAAJI0299R1ZS

	Installation and training at the user's facility, lasting up to three days by trained system engineer Similar equipment should have been installed in government organizations in India during last 2-3 years.	
8	Accessories/Optional: (i) Cuvette pairs 3.5 mL- 2 Nos (ii) Solid Sample Accessory. (iii) Tunable picosecond light source (supercontinuum laser, wavelength range: -360-2000 nm, -200 ps, rep-rate: 10 kHz to 1 MHz)	Accessories/Optional: The following have been removed I. Cuvette pairs 3.5 mL- 2 Nos II. Solid Sample Accessory.
9	Warranty 3-Years	Warranty No change

Anwar Sabath
16/12/15

